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
DENSITY MEASURING SYSTEMS

— Density Gradient Tube Method ◆ Direct Reading —

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For measurement of a large number of samples
With up to 6 large glass tubes

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For measurement of a small number of samples
With up to 4 small glass tubes

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FEATURES

1. Density gradient tube method makes density measurements stable and easy.
2. Standard density floats – the heart of this measurement principle – are backed by Shibayama's 30 years of experience.
3. Backed by maintenance service – re-verification of standard density floats and periodic calibration of the main unit (a calibration certificate is issued).
4. Highly accuracy measurement (Group A floats are calibrated to ± 0.0002 g/cc)
5. Wide range of sample shapes (films, plates, pellets, etc.) and materials (solid form).

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APPLICATIONS

1. Density measurements of high molecular weight materials
2. Density measurements of fibers (natural and man-made)
3. Density measurements of resin products
4. Density measurements of products for quality control

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
The standard density floats in the above density gradient tube are shown at 5 cm intervals for explanation purposes. In an actual application, the floats should be located at 10 to 20 cm intervals.

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ABOUT DENSITY GRADIENT LIQUID

A density gradient liquid is composed of continuous layers of liquids of different densities. The lighter density liquids rise to the top layer of the density gradient liquid column, and the heavier density liquids sink to the bottom. the density gradient of the layers between the top and bottom is continuous.

Installation and Use of Density Gradient Tubes

density gradient tubes filled with the liquid is installed in a constant temperature oven in order to keep the liquid column unchanged over a long period of time. Several standard density floats with even density differences are placed in the tube. Read the scale to show the sample position in reference to the nearest floats.

This measuring principle is based on JIS K 7112-1980 Methods for Determining the Density and Specific Gravity of Plastics.

Preparation of Density Gradient Liquid

Density gradient liquid can be prepared by mixing light and heavy density liquids which have affinity.

Place heavy liquid at the bottom of the glass. Place layers of mixture with lighter density over the heavier layers until the density difference between the light liquid is canceled.

This completes the liquid column with desired density gradient.

WHAT ARE STANDARD DENSITY FLOATS?

The standard density floats are glass balls 5 mm in diameter. The density is adjusted to ± 0.0002 g/cc (Group A) of nominal density.

When these floats are placed in a cylinder of density gradient solution, the graduated scales of the cylinder corresponding with the position of the floats equals the known densities of the floats.

Each float is supplied with a verification certificate showing the test instruments carried out according to JIS K 0061-1992 Test Methods for Density and Relative Density of Chemical Products. For details of test, see 'Traceability for Standard density Floats' separately available.

Specifications and Accuracy Groups

Material: Hard glass

Size: Approx. 5 mm in diameter

Accuracy: Group A (± 0.0002 g/cc) Group B (± 0.001 g/cc)

0.8000 to 1.1999 g/cc 0.800 to 1.199 g/cc

1.2000 to 1.5999 g/cc 1.200 to 1.599 g/cc

1.6000 to 2.2000 g/cc 1.600 to 2.200 g/cc

Usage

Several pieces of standard density floats with even density differences are placed in a density gradient solution filled glass cylinder. Place a few samples whose densities are unknown in the cylinder. Read the graduated scales where the samples rest.

Correlation between reading of the scales on the cylinder in relation to the floats and the resting positions of the samples functions as a measure of the density gradient to determine the density of the samples.

- 1) Standard density floats can be manufactured to provide density values of even intervals (for Group A floats, in the range from 0.8000 to 2.2000 g/cc)
- 2) Each standard density float is supplied with a verification certificate with traceability data.
- 3) If resultant readings are entered into a personal computer (with Lotus or Excel software), the density values can be easily calculated.
- 4) A number of samples can be simultaneously measured.

MEASURING PROCEDURES

Place floats at even intervals in the density gradient tube. → Place samples to measure in the tube. → Read the locations of samples with the graduated scale on the tube. → Enter data values to accessory data mathematical calculation software (which operates on Lotus or Excel) → Density values are calculated according to the method of least squares.

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
SOFTWARE

— for Density Gradient Liquid Preparation, Density Calculation and Graph Printout —
Density Calculation

Calculate density values by entering the tube scale readings of the samples boxes in the table.

The least squares method is used for calculation. To apply the least squares method, fulfill the table with necessary values.

- Input the standard density floats tube scale readings and actual density of the standard density floats.

	1	2	3	4
n	1	1	1	1
Scale Reading : X	46.7	42	34.9	29.7
Density : Y	1.3801	1.3848	1.3902	1.3949
Theoretical Value Y	1.380452	1.384376	1.390305	1.394647
	5	6	7	8
	1	1	1	0
	23	17.2	11.2	0
	1.4	1.4051	1.4101	0
	1.400241	1.405084	1.410094	1.419446

Scale Reading: x = Tube scale reading of standard density float

Density Value: y = Actual density of standard density float

When graphing the resultant data, 'x' and 'y' become X-axis (horizontal) and Y-axis (vertical).

'n' is the number of standard density floats. Generally, one float is used for each density value, n = 1.

The Y theoretical value 'y' is obtained by the following table and expressed by the equation $y = ax + b$.

- Using the data collected in 1., the following table can be obtained.

No.	n	x	y	xy	x^2
1	1	46.7	1.3801	64.45067	2180.89
2	1	42	1.3848	58.1616	1764
3	1	34.9	1.3902	48.51798	1218.01
4	1	29.7	1.3949	41.42853	882.09
5	1	23	1.4	32.2	529
6	1	17.2	1.4051	24.16772	295.84
7	1	11.2	1.4101	15.79312	125.44
8	0	0	0	0	0
計	n = 7	$\Sigma x = 204.7$	$\Sigma y = 9.7652$	$\Sigma xy = 284.7196$	$\Sigma x^2 = 6995.27$

By using the values in the total row, the least squares method equation can be obtained.

Equation: $y = ax + b$

$a = -0.00083$ $b = 1.419446$

After obtaining values for a and b, the theoretical value for y can be calculated.

Verify the calculation equation by substituting an arbitrary number for the tube scale reading.

Tube Scale Reading: 42

Density Value (y theoretical value): 1.3844

- Input the tube scale reading for the sample to obtain, a density value can be obtained.

Tube scale reading for the sample

10	20	30	40
----	----	----	----

Density calculated from tube scale reading

1.4111	1.4027	1.3944	1.386
--------	--------	--------	-------

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
Automatic Constant Speed Sample Lift

This unit is used to minimize the disturbance of density gradient liquid column in the tube by automatically lifting a sample up or down at a constant speed. This unit can maintain the effective density gradient two to three times longer than when using the standard accessory float lifting rod. Prolonging the life of the density gradient saves time and measurement costs.

This unit is mounted on the stirring motor base located at the top of the constant temperature oven.

Main Specifications

- Lifting Up/Down Speeds: 0.9 to 7.5 cm/min.
- Control Method: Integrated timer
- Dimensions and Weight: 100 (W) × 200(D) × 335(H) mm, 3.5 kg
- Power Requirements: 100 VAC, 50/60 Hz. 6 W

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Automatic Density Gradient Liquid Maker

Composed of a pump and a controller, this unit is used to make density gradient liquid with a good linearity. Gradient can be made liquid with significantly improved ease and stability compared with the standard accessory manual liquid maker kit which uses a siphon.

This is recommended for less experienced users. This unit is installed on the top of the constant temperature oven.

Main Specifications

- Control Method: Volume sensor
- Tubing Pump: 4-gang type
- Output Volume: 32 to 450 mL/h × 4
- Magnetic Stirrer: 100 to 1500 rpm
- Power Requirements: 100 VAC, 50/60 Hz. 20 W

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Standard Density Floats

- Other than those supplied as standard accessories -

Built-in Type Cooling Circulator

This unit is installed under the constant temperature oven to keep it at those below room temperature.

Main Specifications

- Control Power: 510/540 kcal/h
- Surrounding Temperature Requirements: 5 to 40 °C
- Cooling Water Temperature Range: 5 to 25 °C
- Power Requirements: 100 VAC, 50/60 Hz. 370 W

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ORDERING INFORMATION

1. Measuring temperature
2. Measuring range
3. Range, quantity and accuracy for standard density floats to be supplied as standard accessory
4. Optional item requirements

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SYSTEM COMPOSITION AND MAIN SPECIFICATIONS

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Measuring Constant Temperature Oven

Temperature Range (Type A and B):
 Room temperature +5 to 40 °C
 When optional built-in type cooling circulator is used,
 15 to 40 °C

Temperature Accuracy (Types A and B): ±0.005 °C
 Temperature Control (Types A and B): PID control
 with auto tuning function

Number of Cylinders:	6	4
Capacity:	134 liters	75 liters
Size (W×D×H mm):	580×480×1590	580×480×1260
Weight:	95 kg	75 kg

Power Requirements (Types A and B): 100 VAC 12 A
 When optional units are used, 17 A

Measuring Glass Tubes

Size (L×Inside Dia.):	1100mm×45 φ	730 mm×32 φ
Scale Range:	850 mm	500 mm

(1 mm intervals)

Standard Density Floats

Standard Accessories (0.800 to 1.199 g/cc):
 30 pcs. 10 pcs.

Note: If floats of other densities and accuracy are needed, place an order as optional accessories,

Soft ware (Types A and B)

Medium: 3.5" FD Operated on Lotus or Excel
 1) Density gradient liquid preparation
 2) Density calculation and graph printout

Manual Density Gradient Liquid Maker

consisting of:

Adjusting bottle	2 pcs. (2 L)	2 pcs. (1 L)
Silicone plug	3 pcs.	2 pcs.
Teflon® tube	2 m	2 m
Tygon® tube	2 m	2 m
Roller clamp	2 pcs.	2 pcs.
Safety pipette	1 pc.	1 pc.
Cylinder weight	1 pc.	1 pc.

Standard Accessories and Documents (Types A and B)

- 1) Thermometer (with calibration certificate) ... 1 pc.
- 2) Float lifting rod 1 pc.
- 3) Operation manual 1 copy
- 4) Test data sheet (with calibration certificate) . . 1 copy
- 5) Standard density float test certificate 1 copy
- 6) Warranty (attached to item 4) 1 copy

Recommended Maintenance Service Intervals

1. Re-verification of standard density floats . . 2-year
2. Thermometer 1-year
3. Measuring constant temperature oven . . . 2-year

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SHIBAYAMA SCIENTIFIC CO., LTD.