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
Polymer Film X and Y-Axis Stretching System

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
Stretching mechanism holding a stretching sample film

COMPOSITION

1. Controller

This unit is designed to control the temperature in the air constant temperature chamber, stretching mode, stretching directions and speeds.

2. Main Unit


1) Air Constant Temperature Chamber/Base

This chamber is composed of a lid, temperature insulated housing and mechanism. The housing contains heaters, stirring fan, sample chucks and stretching mechanism. The temperature can be set up to 300 °C as standard. The lid has a window on its front and rear. The base is an iron plate to mount the system except for the controller.

2) Stretching Mechanism (See photos below.)

The mechanism is provided in the chamber, and is composed of two ball screws, and X and Y-axis sample chucks set on the ball screws. The driving motors and gearheads are provided outside the chamber.

The ball screw mechanisms carry the chucks in the X and Y-axis directions at a set speed to stretch the sample. The motor/gearhead units are controlled by the signals from the controller.

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Stretching mechanism set with an optional stretched sample take-out attachment

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3) Stretching Rate Multiplier Section

The revolution of the ball screw is transferred to a counter via a spur gear to count the number of revolutions. This count is converted into a multiplication signal which determines the rate.

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

This is an X-shaped mechanism designed to fix the positions of the sample chucks, and hold the sample in its stretched status.

Note:
This attachment limits the stretching length to 150 mm, max.

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側面図 Side View

2) 60/90 mm Effective Width Sample Chucks

A pair of 30 mm effective width sample chucks is supplied as standard accessory for a net 30 mm square sample. If a stretching of larger samples are needed, optional 60 mm and 90 mm sample chucks are available.


Effective Chuck Width		Max. Stretching Rate (Length)
30 mm	Std. Accessory	6 times (180 mm)
60 mm	Optional Accessory	3 times (180 mm)
90 mm		2 times (180 mm)

3) Stretched Sample Take-out Attachment

3) Stretched Sample Take-out Attachment

This is balancing mechanism to make opening and closing of the heavy chamber lid (approx. 60 kg) easier, smoother and safer.

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1) 400°C Constant Temperature Chamber

This provides higher temperature of up to 100°C over the standard maximum chamber temperature of 300°C. A partition is provided to increase the thermal efficiency.

2) Max. 300 mm Stretching Mechanism

This provides a larger stretching of 300 mm in place of the standard 180 mm. As this special specification will need a larger base area of 1000×900 mm, specify it when ordering.

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1. Controller

Temperature Controller: Digital controller
Temperature Setting: Key switch operation, with a memory function.
Control Switches: For heaters and stirring fan
Stretching Modes: Single-axis and dual simultaneous dual-axes, and individual dual-axes
Control Switches: Heaters and air stirring fan
Safety Devices: Leak current circuit breaker, 20 A

2. Main Unit

1) Air Constant Temperature Chamber

Temperature Range: Room temperature to 300 °C
Temperature Accuracy: Within ±2%
Heaters: Nichrome wire 1 kW and 500 W
Air Stirring Fan: Induction motor driven sirocco fan
Chamber Construction: Two-split type upper/lower shells. The upper shell is a lid which can be swung open on hinges.
Over-temperature Prevention Device: Thermostat

2) Stretching Mechanism

Stretching Method: Ball screw mechanism, made of SUS, M22-P2.5, with molybdenum disulfide coating. Stretching and returning are made by reversing the motor revolution.
Stretching Mechanism Driving Source: 70 W reversible motors and reduction gears for 15 rpm and 100 rpm
Sample Chucks: Effective width: 30 mm, SUS made
One pair of chucks on each of X and Y-axis.
Stretching Speeds: 15, 20, 30, 100, 132 and 200 mm/min
Stretching Rate: Up to 30 mm × 6 times
Size of Samples: 50mm* square with 1 μ to 1 mm thickness
*Net 30 mm plus a 10 mm margin on either end.

3) Stretching Rate Setting Section

Mechanical Counter: 4-digit display. Stretching rate is controlled by the number of ball screw revolutions.
Two counters are provided for each of X and Y-axes for control of stretching and returning length.

3. Power Source/Dimensions/Weight

Power Source: AC 100 V, 50Hz, or, 60 Hz, 17 A

Dimensions and Weight:

Controller: 430 × 430 × 220 (mm), 15 kg

Main Unit: 890 × 790 × 410 (mm), 150kg

We reserve the right to change specifications without notice.

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Main Products

Density Measuring Systems/Automatic Viscometers/
Polymer Film X and Y Axis Stretching System

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SHIBAYAMA SCIENTIFIC CO., LTD.