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### PC-controlled automatic viscometer

Use a PC to control the viscometer. The SS-600-L1 allows continuous measurement: Washing → Sample Injection → Measurement → Data Processing

### Data management by PC

Data storage, searching, and editing can be performed.

### Process monitoring on the PC

The PC-controlled measuring process can be monitored on the PC.

### Measurement down to 1/100 seconds

Stable data can be obtained down to 1/100 seconds with a 1MHz crystal oscillator.

### Serial measurement of up to 20 samples

Up to 20 samples can be measured under the specified measuring conditions overnight.

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### For measurement of a small number of samples

SS-500-L1 is equivalent to SS-600-L1 without the  
automatic sample feeder.

The automatic sample feeder can be additionally  
installde later.

SS-500-L1 and 600-L1 have a common design  
for easy installation of up to two channels.

For more details, descriptive leaflets are available upon request.

The data calculation software for SS-500-L1 and  
SS-600-L1 uses one of the calculation equations  
shown below.

Select one for your application.

We reserve the right to change the equations without notice.

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### Microprocessor-controlled Automatic Viscometer — RS-232C and Data Calculation Software Included —

Based on our 30 years of experience in automatic viscometers Shibayama offers the following optional accessories to expand the capability of your SS-200-STD.

- a) A P' detector for high viscosity measurement.
- b) A high/low viscosity measurement selector for the P' detector equipped system. This allows fast and easy changing of the system setup for different samples.
- c) A tolerance judgment function to improve measurement efficiency and accuracy.
- d) High temperature measurement up to 135°C.
- e) Automatic sample feeder for continuous measurement of many samples over a long period of time.
- f) An additional channel for dual channel simultaneous measurement.
- g) scanning measurement of different viscosity samples with different types of solvent.

Final specifications for your viscometer will be determined after consultation and determination of desired system configuration.

The data calculation software can provide data based on the resultant mean flow-through time. The calculation equation to be applied can be selected from the six equations on page 2.

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## Microprocessor-controlled Automatic Viscometer

The SS-402-L1 is composed of the SS-200-STD plus a dedicated sequencer and a diluter. Although a raw sample liquid injection is manually performed, the rest of the processes – washing, dilution, measurement and drainage – are automatic. For a washing agent, either a solvent or washing solution can be used. Highly accurate dilution (within 0.1 cc) can be obtained with a gas-tight syringe. Although the number of diluting operations is set to 2 as standard, more diluting operations are available upon request. The number of measurements can be selected from 3 to 5. The number of channels can be increased to two.

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
### Main Specifications

	SS-600-L1	SS-200-STD
Measuring Temperature Range	10 to 40℃	
Measuring Time Range	0.01 to 360 sec.	
Measuring Viscosity perature Range	Within $\eta_{rel}$ 2.3	
Control Method	PC	Microprocessor
Number of Measuring Samples	20	1

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## Review of Shibayama Automatic Viscometer Systems

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For functions not shown in the above examples, please contact us.

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We reserve the right to change specifications without notice.

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**SHIBAYAMA SCIENTIFIC CO., LTD.**